


# Understanding the Art of Stenting

## **Easyflex** *Flexible Coronary Stent System*




One of the  
most  
flexible stent  
designs

# Easyflex

Flexible Coronary Stent System

Two opposite articulations between every stent segment guarantee multidirectional flexibility



- Easy side branch access
- Safe and accurate delivery
- 316 LVM stainless steel
- Minimal foreshortening
- Unique star segments
- Low crimping profile
- Biocompatible
- Balloon expandable
- High flexibility
- Good trackability



## Balloon Mounted Stent

Catalog Number	Balloon Diameter	Stent Length
EF 0725	2,50 mm	07 mm
EF 0727	2,75 mm	07 mm
EF 0730	3,00 mm	07 mm
EF 0735	3,50 mm	07 mm
EF 0740	4,00 mm	07 mm
EF 1025	2,50 mm	10 mm
EF 1027	2,75 mm	10 mm
EF 1030	3,00 mm	10 mm
EF 1035	3,50 mm	10 mm
EF 1040	4,00 mm	10 mm
EF 1325	2,50 mm	13 mm
EF 1327	2,75 mm	13 mm
EF 1330	3,00 mm	13 mm
EF 1335	3,50 mm	13 mm
EF 1340	4,00 mm	13 mm
EF 1725	2,50 mm	17 mm
EF 1727	2,75 mm	17 mm
EF 1730	3,00 mm	17 mm
EF 1735	3,50 mm	17 mm
EF 1740	4,00 mm	17 mm
EF 2025	2,50 mm	20 mm
EF 2027	2,75 mm	20 mm
EF 2030	3,00 mm	20 mm
EF 2035	3,50 mm	20 mm
EF 2040	4,00 mm	20 mm
EF 2425	2,50 mm	24 mm
EF 2427	2,75 mm	24 mm
EF 2430	3,00 mm	24 mm
EF 2435	3,50 mm	24 mm
EF 2440	4,00 mm	24 mm
EF 2825	2,50 mm	28 mm
EF 2827	2,75 mm	28 mm
EF 2830	3,00 mm	28 mm
EF 2835	3,50 mm	28 mm
EF 2840	4,00 mm	28 mm
EF 3325	2,50 mm	33 mm
EF 3327	2,75 mm	33 mm
EF 3330	3,00 mm	33 mm
EF 3335	3,50 mm	33 mm
EF 3340	4,00 mm	33 mm
EF 3825	2,50 mm	38 mm
EF 3827	2,75 mm	38 mm
EF 3830	3,00 mm	38 mm
EF 3835	3,50 mm	38 mm
EF 3840	4,00 mm	38 mm

## Bare Stent

Catalog Number	Stent Length
EF 0700B	07 mm
EF 1000B	10 mm
EF 1300B	13 mm
EF 1700B	17 mm
EF 2000B	20 mm
EF 2400B	24 mm
EF 2800B	28 mm
EF 3300B	33 mm
EF 3800B	38 mm

distributed by

**iramed** medical devices

iramed GmbH, Germany  
Im End 39, D-41844 Wegberg-Arsbeck  
Tel. +49-2436-588, Fax +49-2436-1091  
E-mail: iramed@aol.com, www.iramed.com