




Understanding the Art of Stenting

MovyTax DLC Paclitaxel Eluting Stent System

- **Coated by Biosteel® with DLC (Diamond Like Carbon)** for longterm biocompatibility and covered with Polymer & Paclitaxel for slow release within 30 days
- **Star shaped segments with flexible links** to follow vessel movements



Taxus brevifolia



MovyTax

DLC Paclitaxel Eluting Stent System

Balloon Mounted Stent

Catalog Number	Balloon Diameter	Stent Length
MT 0725	2,50 mm	07 mm
MT 0727	2,75 mm	07 mm
MT 0730	3,00 mm	07 mm
MT 0735	3,50 mm	07 mm
MT 0740	4,00 mm	07 mm
MT 1025	2,50 mm	10 mm
MT 1027	2,75 mm	10 mm
MT 1030	3,00 mm	10 mm
MT 1035	3,50 mm	10 mm
MT 1040	4,00 mm	10 mm
MT 1325	2,50 mm	13 mm
MT 1327	2,75 mm	13 mm
MT 1330	3,00 mm	13 mm
MT 1335	3,50 mm	13 mm
MT 1340	4,00 mm	13 mm
MT 1725	2,50 mm	17 mm
MT 1727	2,75 mm	17 mm
MT 1730	3,00 mm	17 mm
MT 1735	3,50 mm	17 mm
MT 1740	4,00 mm	17 mm
MT 2025	2,50 mm	20 mm
MT 2027	2,75 mm	20 mm
MT 2030	3,00 mm	20 mm
MT 2035	3,50 mm	20 mm
MT 2040	4,00 mm	20 mm
MT 2425	2,50 mm	24 mm
MT 2427	2,75 mm	24 mm
MT 2430	3,00 mm	24 mm
MT 2435	3,50 mm	24 mm
MT 2440	4,00 mm	24 mm
MT 2825	2,50 mm	28 mm
MT 2827	2,75 mm	28 mm
MT 2830	3,00 mm	28 mm
MT 2835	3,50 mm	28 mm
MT 2840	4,00 mm	28 mm
MT 3325	2,50 mm	33 mm
MT 3327	2,75 mm	33 mm
MT 3330	3,00 mm	33 mm
MT 3335	3,50 mm	33 mm
MT 3340	4,00 mm	33 mm
MT 3825	2,50 mm	38 mm
MT 3827	2,75 mm	38 mm
MT 3830	3,00 mm	38 mm
MT 3835	3,50 mm	38 mm
MT 3840	4,00 mm	38 mm

MovyTax

DLC Paclitaxel Eluting Stent System



Target dose: 1 µg/mm ²
Biosteel coating technology
Material: Stainless Steel
DLC (diamond like carbon coated) covered with polymer & drug
Therapeutic release within 28 days
High mechanical stability of the coating

Technical specifications

Drug:	Paclitaxel 1µg/mm ²
Polymer:	Inert polymer compound/fully resorbable
Release kinetic:	Slow release within 10-28 days
Coating thickness:	1-2 µm
Stent surface finishing:	DLC PII (plasma ion implantation)

Bare Stent

Catalog Number	Stent Length
MT 0700B	07 mm
MT 1000B	10 mm
MT 1300B	13 mm
MT 1700B	17 mm
MT 2000B	20 mm
MT 2400B	24 mm
MT 2800B	28 mm
MT 3300B	33 mm
MT 3800B	38 mm



Post stenting 180 days

Strut Composition



Longitudinal cut through stent strut